



CD47 Rabbit anti-Human Polyclonal (aa23-140) Antibody - LS-B10955 - LSBio	
CatalogID:	LS-B10955
Validation:	This antibody replaces catalog number LS-C178715. It has been validated for use in the following assays: IHC-P.
Target:	CD47 molecule
Synonyms:	CD47 Antibody, CD47 antigen Antibody, CD47 molecule Antibody, Integrin associated protein Antibody, Integrin-associated protein Antibody, Leukocyte surface antigen CD47 Antibody, MER6 Antibody, Protein MER6 Antibody, Rh-related antigen Antibody, OA3 Antibody, CD47 glycoprotein Antibody, IAP Antibody
Host	CD47 antibody was produced in Rabbit
Clonality:	Polyclonal
Immunogen Species:	CD47 antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	CD47 antibody was raised against a recombinant human CD47 protein 23-140aa (BC010016).
Specificity:	Human CD47.
Epitope:	aa23-140
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.4, 0.02% sodium azide
Recommended Storage:	Store at -20°C for up to one year.
Uses:	IHC - Paraffin (2.5 - 5 µg/ml), Western blot (1:1000 - 1:2000), ELISA (1:20000 - 1:80000) (Optimal dilution to be determined by the researcher)
Size:	50 µg
Concentration:	1 mg/ml

Immunohistochemistry Image:



Human Prostate: Formalin-Fixed, Paraffin-Embedded (FFPE)

Western Blot Image:



**Western Blot analysis of
Jurkat cell Lysate with anti
-CD47 antibody**

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/23/2014

© 2014 LifeSpan BioSciences