


HDGFRP2 Rabbit anti-Human Polyclonal (aa51-79) Antibody - LS-B10977 - LSBio


CatalogID:	LS-B10977
Validation:	This antibody replaces catalog number LS-C157802. It has been validated for use in the following assays: IHC-P.
Target:	hepatoma-derived growth factor-related protein 2 (HDGFRP2)
Synonyms:	HDGFRP2 Antibody, HDGF-2 Antibody, HDGF2 Antibody, HRP-2 Antibody
Host	HDGFRP2 antibody was produced in Rabbit
Clonality:	Polyclonal
Immunogen Species:	HDGFRP2 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	HDGFRP2 antibody was raised against kLH-conjugated synthetic peptide from N-terminal region of human HDGR2.
Specificity:	Human HDGFRP2
Epitope:	aa51-79
Reactivity:	Human
Purification:	Immunoaffinity purified
Presentation:	PBS, 0.09% sodium azide
Recommended Storage:	Maintain refrigerated at 2-8°C for up to 6 months. For long term storage store at -20°C.
Uses:	IHC - Paraffin (1:100), Western blot (1:1000) (Optimal dilution to be determined by the researcher)
Size:	200 µl

Immunohistochemistry Image:



Human Testis: Formalin-Fixed, Paraffin-Embedded (FFPE)

Western Blot Image:



HDGR2 Antibody (N-term) Western blot of HepG2 cell line lysates (35 ug/lane). This demonstrates the HDGR2 antibody detected the HDGR2 protein (arrow).

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/23/2014

© 2014 LifeSpan BioSciences