


SALL4 Rabbit anti-Mouse Polyclonal Antibody - LS-B138 - LSBio

CatalogID:	LS-B138
Validation:	This antibody replaces catalog number LS-C3306. It has been validated for use in the following assays: IHC.
Target:	spalt-like transcription factor 4 (SALL4)
Synonyms:	SALL4 Antibody, DRRS Antibody, DJ1112F19.1 Antibody, HSAL4 Antibody, Sal-like 4 (Drosophila) Antibody, Sal-like protein 4 Antibody, Zinc finger protein SALL4 Antibody, ZNF797 Antibody, Sal (Drosophila)-like 4 Antibody, Zinc finger protein 797 Antibody
Host	SALL4 antibody was produced in Rabbit
Clonality:	Polyclonal
Immunogen Species:	SALL4 antibody was raised against Mouse
Antigen Type:	Synthetic peptide
Immunogen:	SALL4 antibody was raised against synthetic peptide from mouse SALL4.
Reactivity:	Mouse, Human
Purification:	Protein G purified
Presentation:	PBS, 0.09% sodium azide.
Recommended Storage:	Long term: -20°C; Short term: +4°C. Avoid repeat freeze-thaw cycles.
Usage Summary:	Immunohistochemistry: LS-B138 was validated for use in immunohistochemistry on a panel of 21 formalin-fixed, paraffin-embedded (FFPE) human tissues after heat induced antigen retrieval in pH 6.0 citrate buffer. After incubation with the primary antibody, slides were incubated with biotinylated secondary antibody, followed by alkaline phosphatase-streptavidin and chromogen. The stained slides were evaluated by a pathologist to confirm staining specificity. The optimal working concentration for LS-B138 was determined to be 5-10 ug/ml.
Uses:	IHC - Paraffin (5 - 10 µg/ml), ELISA (1:000 - 1:1000) (Optimal dilution to be determined by the researcher)
Size:	50 µg
Concentration:	1 mg/ml

Immunohistochemistry Image:



Anti-SALL4 antibody IHC of human skeletal muscle. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval. Antibody LS-B138 concentration 10 ug/ml.

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/23/2014

© 2014 LifeSpan BioSciences