



TNIP2 / ABIN-2 Rabbit anti-Human Polyclonal Antibody - LS-B2082 - LSBio	
CatalogID:	LS-B2082
Validation:	This antibody replaces catalog number LS-C613. It has been validated for use in the following assays: IHC-P.
Target:	TNFAIP3 interacting protein 2 (TNIP2)
Synonyms:	TNIP2 Antibody, ABIN-2 Antibody, ABIN2 Antibody, FLIP1 Antibody, KLIP Antibody, TNFAIP3-interacting protein 2 Antibody, TNFAIP3 interacting protein 2 Antibody
Host	TNIP2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Immunogen Species:	TNIP2 / ABIN-2 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	TNIP2 / ABIN-2 antibody was raised against synthetic peptide from human TNIP2.
Specificity:	Synthetic peptides of human ABIN-2. It cross-reacts with mouse ABIN-2.
Reactivity:	Human, Mouse
Purification:	Protein G purified
Presentation:	PBS, 0.05% BSA, 0.05% sodium azide
Recommended Storage:	Short term 4°C, long term aliquot and store at -20°C, avoid freeze thaw cycles.
Usage Summary:	Immunohistochemistry: LS-B2082 was validated for use in immunohistochemistry on a panel of 21 formalin-fixed, paraffin-embedded (FFPE) human tissues after heat induced antigen retrieval in pH 6.0 citrate buffer. After incubation with the primary antibody, slides were incubated with biotinylated secondary antibody, followed by alkaline phosphatase-streptavidin and chromogen. The stained slides were evaluated by a pathologist to confirm staining specificity. The optimal working concentration for LS-B2082 was determined to be 5 ug/ml.
Uses:	IHC - Paraffin (5 µg/ml), Western blot (0.5 - 2 µg/ml) (Optimal dilution to be determined by the researcher)
Size:	50 µg
Concentration:	0.5 mg/ml

Immunohistochemistry Image:



Anti-TNIP2 antibody IHC of human brain cortex. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval. Antibody LS-B2082 concentration 5 ug/ml.

Western Blot Image:



Western blot of ABIN-2 in A) human brain, B) mouse brain and C) rat brain lysate using antibody at 3 ug/ml.

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/23/2014

© 2014 LifeSpan BioSciences