


MMP9 / Gelatinase B Mouse anti-Human Monoclonal Antibody - LS-B2371 - LSBio	
CatalogID:	LS-B2371
Validation:	This antibody replaces catalog number LS-C8809. It has been validated for use in the following assays: IHC-P.
Target:	matrix metalloproteinase 9 (gelatinase B, 92kDa gelatinase, 92kDa type IV collagenase) (MMP9)
Synonyms:	MMP9 Antibody, 92kd gelatinase Antibody, 92 kd type iv collagenase Antibody, 92 kDa gelatinase Antibody, 92 kDa type IV collagenase Antibody, Matrix metalloproteinase-9 Antibody, MMP-9 Antibody, Type V collagenase Antibody, CLG4B Antibody, Gelatinase B Antibody, GELB Antibody, Macrophage gelatinase Antibody, MANDP2 Antibody
Family / Subfamily:	Protease / Metalloproteinase M10A
Host	MMP9 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Immunogen Species:	MMP9 / Gelatinase B antibody was raised against Human
Antigen Type:	Purified protein
Immunogen:	MMP9 / Gelatinase B antibody was raised against purified human gelatinase B. Cellular Localization: Peripheral cytoplasmic.
Specificity:	Recognizes human MMP-9 at 92kD and 86kD (non-reducing conditions) which are identified as pro (latent) and active forms of matrix metalloproteinase-9 (MMP-9; also known as 92kD and collagenase IV or gelatinase B). Reacts with the complex (~200kD) of MMP-9 and tissue inhibitor of metalloproteinases-1 (TIMP-1)1. Does not cross-react with pro and active forms of MMP-2.
Reactivity:	Human
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 0.2% BSA, 0.09% sodium azide.
Recommended Storage:	Long term: -20°C; Short term: +4°C. Avoid repeat freeze-thaw cycles.
Usage Summary:	Immunohistochemistry: LS-B2371 was validated for use in immunohistochemistry on a panel of 21 formalin-fixed, paraffin-embedded (FFPE) human tissues after heat induced antigen retrieval in pH 6.0 citrate buffer. After incubation with the primary antibody, slides were incubated with biotinylated secondary antibody, followed by alkaline phosphatase-streptavidin and chromogen. The stained slides were evaluated by a pathologist to confirm staining specificity. The optimal working concentration for LS-B2371 was determined to be 5 ug/ml. Western Blot: 1 ug/ml for 2hrs at RT under non-reducing conditions. Immunoprecipitation: 2 ug/mg protein lysate. Use Protein G. Inhibition: 2-4 ug/ml.
Uses:	IHC - Paraffin (5 ug/ml), IHC - Frozen (1:100), Western blot, Immunoprecipitation (Optimal dilution to be determined by the researcher)
Size:	50 ug
Concentration:	0.2 mg/ml

Immunohistochemistry Image:



Anti-MMP-9 antibody IHC of human spleen, macrophages. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval. Antibody LS-B2371 concentration 5 ug/ml.

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/23/2014

© 2014 LifeSpan BioSciences