


PAWR / Par-4 Rabbit anti-Human Polyclonal (N-Terminus) Antibody - LS-B267 - LSBio

CatalogID:	LS-B267
Validation:	This antibody replaces catalog number LS-C2825. It has been validated for use in the following assays: IHC.
Target:	PRKC, apoptosis, WT1, regulator (PAWR)
Synonyms:	PAWR Antibody, PAR4 Antibody, WT1-interacting protein Antibody, Par-4 Antibody, Transcriptional repressor PAR4 Antibody
Host	PAWR antibody was produced in Rabbit
Clonality:	Polyclonal
Immunogen Species:	PAWR / Par-4 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	PAWR / Par-4 antibody was raised against synthetic peptides corresponding to the N-terminal sequence (MATGGYRTSSGLGGS) and C-terminal sequence (QENKTLKVVGQLTR) of human PAR4. Percent identity by BLAST analysis: Human, Gorilla, Gibbon, Monkey (100%); Panda (93%); Rabbit (87%); Opossum (80%).
Specificity:	PAR4.
Epitope:	N-Terminus
Reactivity:	Human, Gorilla, Gibbon, Monkey
Purification:	Affinity purified
Presentation:	PBS, 0.02% sodium azide.
Recommended Storage:	Long term: -20°C; Short term: +4°C. Avoid repeat freeze-thaw cycles.
Usage Summary:	Immunohistochemistry: LS-B267 was validated for use in immunohistochemistry on a panel of 21 formalin-fixed, paraffin-embedded (FFPE) human tissues after heat induced antigen retrieval in pH 6.0 citrate buffer. After incubation with the primary antibody, slides were incubated with biotinylated secondary antibody, followed by alkaline phosphatase-streptavidin and chromogen. The stained slides were evaluated by a pathologist to confirm staining specificity. The optimal working concentration for LS-B267 was determined to be 10 ng/ml.
Uses:	IHC - Paraffin (10 µg/ml), Western blot (1:500 - 1:1000) (Optimal dilution to be determined by the researcher)
Size:	50 µl
Concentration:	1 mg/ml

Immunohistochemistry Image:



Anti-PAWR / PAR4 antibody IHC of human lung. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval. Antibody LS-B267 concentration 10 ug/ml.

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/23/2014

© 2014 LifeSpan BioSciences