


HIF3A / HIF3-Alpha Rabbit anti-Human Polyclonal (N-Terminus) Antibody - LS-B287 - LSBio	
<b>CatalogID:</b>	LS-B287
<b>Validation:</b>	This antibody replaces catalog number LS-C2566. It has been validated for use in the following assays: IHC.
<b>Target:</b>	hypoxia inducible factor 3, alpha subunit (HIF3A)
<b>Synonyms:</b>	HIF3A Antibody, BHLHe17 Antibody, HIF-3A Antibody, IPAS Antibody, HIF3-alpha-1 Antibody, Inhibitory PAS domain protein Antibody, PASD7 Antibody, HIF-3-alpha Antibody, HIF3-alpha Antibody, Member of PAS protein 7 Antibody, MOP7 Antibody
<b>Host</b>	HIF3A antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Immunogen Species:</b>	HIF3A / HIF3-Alpha antibody was raised against Human
<b>Immunogen:</b>	HIF3A / HIF3-Alpha antibody was raised against bacterially expressed human HIF-3 alpha (N-terminus)
<b>Specificity:</b>	HIF-3 alpha.
<b>Epitope:</b>	N-Terminus
<b>Reactivity:</b>	Human, Mouse
<b>Purification:</b>	Immunoaffinity purified
<b>Presentation:</b>	Tris-citrate/phosphate buffer, pH 7-8, 0.1% sodium azide.
<b>Recommended Storage:</b>	+4°C, avoid freezing
<b>Usage Summary:</b>	Immunohistochemistry: LS-B287 was validated for use in immunohistochemistry on a panel of 21 formalin-fixed, paraffin-embedded (FFPE) human tissues after heat induced antigen retrieval in pH 6.0 citrate buffer. After incubation with the primary antibody, slides were incubated with biotinylated secondary antibody, followed by alkaline phosphatase-streptavidin and chromogen. The stained slides were evaluated by a pathologist to confirm staining specificity. The optimal working concentration for LS-B287 was determined to be 5 ug/ml.
<b>Uses:</b>	IHC - Paraffin (5 µg/ml), Western blot (1:100 - 1:200) (Optimal dilution to be determined by the researcher)
<b>Size:</b>	50 µl
<b>Concentration:</b>	1 mg/ml

**Immunohistochemistry Image:**



Anti-HIF3A antibody IHC of human lung. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval. Antibody LS-B287 concentration 5 ug/ml.

**Requested From:**

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/23/2014

© 2014 LifeSpan BioSciences