


TNFAIP3 / A20 Goat anti-Human Polyclonal Antibody - LS-B3849 - LSBio	
CatalogID:	LS-B3849
Validation:	This antibody replaces catalog number LS-C108266. It has been validated for use in the following assays: IHC-Fr.
Target:	tumor necrosis factor, alpha-induced protein 3 (TNFAIP3)
Synonyms:	TNFAIP3 Antibody, A20 Antibody, TNFA1P2 Antibody, TNF alpha-induced protein 3 Antibody, OTUD7C Antibody, Zinc finger protein A20 Antibody
Host	TNFAIP3 antibody was produced in Goat
Clonality:	Polyclonal
Immunogen Species:	TNFAIP3 / A20 antibody was raised against Human
Immunogen:	TNFAIP3 / A20 antibody was raised against synthetic peptide from human TNFAIP3 / A20. Percent identity by BLAST analysis: Human, Gorilla, Gibbon, Monkey, Marmoset, Rabbit (100%); Mouse, Hamster (93%); Bat, Bovine, Pig (86%).
Specificity:	Human TNFAIP3.
Reactivity:	Human, Gorilla, Gibbon, Monkey, Rabbit
Predicted Reactivity:	Mouse, Hamster
Purification:	Immunoaffinity purified
Presentation:	Tris-buffered saline, pH 7.3, 0.5% BSA, 0.02% sodium azide
Recommended Storage:	Store at -20°C. Minimize freezing and thawing.
Usage Summary:	Immunohistochemistry: LS-B3849 was validated for use in immunohistochemistry on a panel of 21 formalin-fixed, paraffin-embedded (FFPE) human tissues after heat induced antigen retrieval in pH 6.0 citrate buffer. After incubation with the primary antibody, slides were incubated with biotinylated secondary antibody, followed by alkaline phosphatase-streptavidin and chromogen. The stained slides were evaluated by a pathologist to confirm staining specificity. The optimal working concentration for LS-B3849 was determined to be 3.75ug/ml.
Uses:	IHC - Paraffin (3.75 µg/ml) (Optimal dilution to be determined by the researcher)
Size:	50 µg
Concentration:	0.5 mg/ml

Immunohistochemistry Image:



Anti-TNFAIP3 antibody IHC of human small intestine. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval. Antibody LS-B3849 concentration 3.75 ug/ml.

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/23/2014

© 2014 LifeSpan BioSciences