


RTN3 / Reticulon 3 Mouse anti-Human Monoclonal (10A8) Antibody - LS-B4116 - LSBio


CatalogID:	LS-B4116
Validation:	This antibody replaces catalog number LS-C108727. It has been validated for use in the following assays: IHC-P.
Target:	reticulon 3 (RTN3)
Synonyms:	RTN3 Antibody, ASYIP Antibody, ASY interacting protein Antibody, HAP Antibody, NSP-like protein II Antibody, NSPL2 Antibody, RTN3-A1 Antibody, NSP-like protein 2 Antibody, Reticulon-3 Antibody, Homolog of ASY protein Antibody, NSPLII Antibody, Reticulon 3 Antibody
Family / Subfamily:	Reticulon
Host	RTN3 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	10A8
Immunogen Species:	RTN3 / Reticulon 3 antibody was raised against Human
Immunogen:	RTN3 / Reticulon 3 antibody was raised against purified recombinant fragment of RTN3 expressed in E. coli.
Specificity:	Human RTN3
Reactivity:	Human
Purification:	Ascites
Presentation:	Ascites fluid, 0.03% sodium azide
Recommended Storage:	Long term: -20°C; Short term: +4°C. Avoid repeat freeze-thaw cycles.
Uses:	IHC - Paraffin (1:200), ELISA (1:10000) (Optimal dilution to be determined by the researcher)
Size:	50 µl

Immunohistochemistry Image:



Anti-RTN3 / Reticulon 3 antibody IHC of human liver. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval. Antibody LS-B4116 dilution 1:200.

Immunohistochemistry Image:



Anti-RTN3 / Reticulon 3 antibody IHC of human kidney. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval. Antibody LS-B4116 dilution 1:200.

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/23/2014

© 2014 LifeSpan BioSciences