


AGXT / SPT Mouse anti-Human Monoclonal (AT2T4) Antibody - LS-B4187 - LSBio

CatalogID:	LS-B4187
Validation:	This antibody replaces catalog number LS-C113041. It has been validated for use in the following assays: IHC-P.
Target:	alanine-glyoxylate aminotransferase (AGXT)
Synonyms:	AGXT Antibody, AGXT1 Antibody, AGT1 Antibody, PH1 Antibody, SPT Antibody, SPAT Antibody, TLH6 Antibody
Host	AGXT antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2b,k
Clone Name:	AT2T4
Immunogen Species:	AGXT / SPT antibody was raised against Human
Reactivity:	Human
Purification:	Protein G purified
Presentation:	Phosphate buffered saline, pH 7.4, 0.1% sodium azide
Recommended Storage:	+4°C or -20°C, Avoid repeated freezing and thawing.
Uses:	IHC - Paraffin (5 µg/ml) (Optimal dilution to be determined by the researcher)
Size:	50 µl
Concentration:	1 mg/ml

Immunohistochemistry Image:



Anti-AGXT antibody IHC of human liver. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval. Antibody LS-B4187 concentration 5 ug/ml.

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/23/2014

© 2014 LifeSpan BioSciences