


SFTPB / Surfactant Protein B Mouse anti-Human Monoclonal (SPB01) Antibody - LS-B5172 - LSBio	
CatalogID:	LS-B5172
Validation:	This antibody replaces catalog number LS-C88382. It has been validated for use in the following assays: IHC-P.
Target:	surfactant protein B (SFTPB)
Synonyms:	SFTPB Antibody, 6 kDa protein Antibody, PSP-B Antibody, SFTP3 Antibody, SP-B Antibody, SFTB3 Antibody, SMDP1 Antibody, Surfactant protein B Antibody
Host	SFTPB antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a,k
Clone Name:	SPB01
Immunogen Species:	SFTPB / Surfactant Protein B antibody was raised against Human
Immunogen:	SFTPB / Surfactant Protein B antibody was raised against recombinant pro-surfactant protein B.
Reactivity:	Human, Rat
Purification:	Protein A purified
Presentation:	10mM PBS, pH 7.4, with 0.2% BSA and 0.09% Sodium Azide
Recommended Storage:	Stable for 24 months when stored at 2-8°C.
Uses:	IHC - Paraffin (20 µg/ml), Immunofluorescence (Optimal dilution to be determined by the researcher)
Size:	50 µg
Concentration:	0.2 mg/ml

Immunohistochemistry Image:



Anti-SFTPB / SP-B antibody IHC of human lung. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval. Antibody LS-B5172 concentration 20 ug/ml.

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences