


**PAX9 Mouse anti-Human Monoclonal (3B8) Antibody - LS-B6405 - LSBio**


<b>CatalogID:</b>	LS-B6405
<b>Validation:</b>	This antibody replaces catalog number LS-C139130. It has been validated for use in the following assays: IHC-P.
<b>Target:</b>	paired box 9 (PAX9)
<b>Synonyms:</b>	PAX9 Antibody, Paired box 9 Antibody, Paired box protein Pax-9 Antibody, STHAG3 Antibody, Paired box gene 9 Antibody, Paired domain gene 9 Antibody
<b>Family / Subfamily:</b>	Transcription factor / not assigned-Transcription factor
<b>Host</b>	PAX9 antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a,k
<b>Clone Name:</b>	3B8
<b>Immunogen Species:</b>	PAX9 antibody was raised against Human
<b>Antigen Type:</b>	Recombinant protein
<b>Immunogen:</b>	PAX9 antibody was raised against recombinant protein
<b>Specificity:</b>	Human PAX9
<b>Reactivity:</b>	Human
<b>Purification:</b>	Protein A purified
<b>Presentation:</b>	PBS, pH 7.2
<b>Recommended Storage:</b>	Short term 4°C, long term aliquot and store at -20°C, avoid freeze thaw cycles.
<b>Usage Summary:</b>	Western Blot: Transfected lysate, Recombinant protein.
<b>Uses:</b>	IHC - Paraffin (5 µg/ml), Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	50 µg
<b>Concentration:</b>	0.5 mg/ml

**Immunohistochemistry Image:**



Anti-PAX9 antibody IHC of human uterus, endometrium. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval. Antibody LS-B6405 concentration 5 ug/ml.

**Western Blot Image:**



Western blot of PAX9 expression in transfected 293T cell line by PAX9 monoclonal antibody, clone 3B8.

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences