


FOXQ1 Mouse anti-Human Monoclonal (3E3) Antibody - LS-B6410 - LSBio


CatalogID:	LS-B6410
Validation:	This antibody replaces catalog number LS-C139140. It has been validated for use in the following assays: IHC-P.
Target:	forkhead box Q1 (FOXQ1)
Synonyms:	FOXQ1 Antibody, Forkhead box Q1 Antibody, HFH-1 Antibody, HNF-3/forkhead-like protein 1 Antibody, HFH1 Antibody, Forkhead box protein Q1 Antibody
Host	FOXQ1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2b,k
Clone Name:	3E3
Immunogen Species:	FOXQ1 antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	FOXQ1 antibody was raised against recombinant protein
Specificity:	Human FOXQ1
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2
Recommended Storage:	Short term 4°C, long term aliquot and store at -20°C, avoid freeze thaw cycles.
Usage Summary:	Sandwich ELISA: Recombinant protein.
Uses:	IHC - Paraffin (5 µg/ml), ELISA (Optimal dilution to be determined by the researcher)
Size:	50 µg

Immunohistochemistry Image:



Anti-FOXQ1 antibody IHC of human small intestine. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval. Antibody LS-B6410 concentration 5 ug/ml.

ELISA Image:



Detection limit for recombinant GST tagged FOXQ1 is 0.3 ng/ml as a capture antibody.

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences