


PHB / Prohibitin Rabbit anti-Human Polyclonal Antibody - LS-B7282 - LSBio

CatalogID:	LS-B7282
Validation:	This antibody replaces catalog number LS-C146959. It has been validated for use in the following assays: IHC-P.
Target:	prohibitin (PHB)
Synonyms:	PHB Antibody, PHB1 Antibody, Prohibitin Antibody
Host	PHB antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Immunogen Species:	PHB / Prohibitin antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	PHB / Prohibitin antibody was raised against synthetic peptide surrounding amino acid 256 of human Prohibitin
Specificity:	Human PHB / Prohibitin
Reactivity:	Human, Mouse, Rat
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 30% glycerol, 0.5% BSA, 0.01% Thimerosal
Recommended Storage:	Long term: -80°C; Short term: -20°C
Uses:	IHC - Paraffin (5 µg/ml), Immunofluorescence, Western blot (0.5 - 4 µg/ml), Immunoprecipitation (10 - 20 µg/ml) (Optimal dilution to be determined by the researcher)
Size:	50 µg
Concentration:	0.5 mg/ml

Immunohistochemistry Image:



Anti-PHB / Prohibitin antibody IHC of human colon. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval. Antibody LS-B7282 concentration 5 ug/ml.

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences