



VEGFA / VEGF Rabbit anti-Rat Polyclonal Antibody - LS-B7747 - LSBio	
CatalogID:	LS-B7747
Validation:	This antibody replaces catalog number LS-C147117. It has been validated for use in the following assays: IHC-P.
Target:	vascular endothelial growth factor A (VEGFA)
Synonyms:	VEGFA Antibody, VPF Antibody, Vascular permeability factor Antibody, VEGF Antibody, VEGF-A Antibody, MVCD1 Antibody
Family / Subfamily:	PDGF / not assigned-PDGF
Host	VEGFA antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Immunogen Species:	VEGFA / VEGF antibody was raised against Rat
Antigen Type:	Recombinant protein
Immunogen:	VEGFA / VEGF antibody was raised against recombinant rat VEGF expressed from E. coli
Specificity:	Rat VEGFA / VEGF
Reactivity:	Rat, Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 30% glycerol, 0.5% BSA, 0.01% Thimerosal
Recommended Storage:	Short term 4°C, long term aliquot and store at -20°C, avoid freeze thaw cycles.
Uses:	IHC - Paraffin (5 µg/ml), Western blot (0.5 - 4 µg/ml) (Optimal dilution to be determined by the researcher)
Size:	50 µg
Concentration:	0.5 mg/ml

Immunohistochemistry Image:



Anti-VEGF antibody IHC of human lung, vascular endothelium. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval. Antibody LS-B7747 dilution 5 ug/ml.

Western Blot Image:



Western blot of VEGF antibody LS-B7747.

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences