


A1BG / GAB Mouse anti-Human Monoclonal Antibody - LS-B7846 - LSBio


CatalogID:	LS-B7846
Validation:	This antibody replaces catalog number LS-C62080. It has been validated for use in the following assays: IHC-P.
Target:	alpha-1-B glycoprotein (A1BG)
Synonyms:	A1BG Antibody, A1B Antibody, ABG Antibody, GAB Antibody, HYST2477 Antibody, Alpha-1-B glycoprotein Antibody, Alpha-1B-glycoprotein Antibody
Host	A1BG antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1,k
Immunogen Species:	A1BG / GAB antibody was raised against Human
Immunogen:	A1BG / GAB antibody was raised against protein purified from Human plasma.
Reactivity:	Human
Purification:	Ammonium sulfate precipitation
Presentation:	HEPES with 0.15M NaCl, 0.01% BSA, 0.03% sodium azide, and 50% glycerol
Recommended Storage:	Long term: -20°C; Short term: -20°C
Uses:	IHC - Paraffin (1:200), Western blot (1:2000), ELISA (Optimal dilution to be determined by the researcher)
Size:	50 µl

Immunohistochemistry Image:



Anti-A1BG / GAB antibody IHC of human colon, submucosal vessels. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval. Antibody LS-B7846 dilution 1:200.

Western Blot Image:



Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/24/2014	
© 2014 LifeSpan BioSciences	