


Hemoglobin Mouse anti-Human Monoclonal Antibody - LS-B7868 - LSBio

CatalogID:	LS-B7868
Validation:	This antibody replaces catalog number LS-C83309. It has been validated for use in the following assays: IHC-P.
Target:	Hemoglobin
Host	Hemoglobin antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Immunogen Species:	Hemoglobin antibody was raised against Human
Antigen Type:	Purified protein
Immunogen:	Hemoglobin antibody was raised against purified human hemoglobin.
Specificity:	Does not cross-react with hemoglobin of the following species: bovine, chicken, equine, ovine, or porcine.
Reactivity:	Human
Purification:	Protein G purified
Presentation:	PBS, pH 7.4
Recommended Storage:	Short term 4°C, long term aliquot and store at -20°C, avoid freeze thaw cycles.
Usage Summary:	May be used in immunoassays to detect and quantitate human hemoglobin.
Uses:	IHC - Paraffin (10 µg/ml), ELISA (Optimal dilution to be determined by the researcher)
Size:	250 µg
Concentration:	1.55 mg/ml

Immunohistochemistry Image:



Anti-Hemoglobin antibody IHC of human placenta, erythrocytes. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval. Antibody LS-B7868 dilution 10 ug/ml.

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences