


EIF3A / EIF3 Rabbit anti-Human Polyclonal (aa550-600) Antibody - LS-B7957 - LSBio


CatalogID:	LS-B7957
Validation:	This antibody replaces catalog number LS-C148506. It has been validated for use in the following assays: IHC-P.
Target:	eukaryotic translation initiation factor 3, subunit A (EIF3A)
Synonyms:	EIF3A Antibody, Centrosomin homolog Antibody, Cytoplasmic protein p167 Antibody, EIF3S10 Antibody, EIF3 p167 Antibody, EIF3 p180 Antibody, EIF-3-theta Antibody, EIF3 p185 Antibody, EIF3, p180 subunit Antibody, KIAA0139 Antibody, p167 Antibody, p180 Antibody, p185 Antibody, TIF32 Antibody, EIF3 Antibody, EIF3-p170 Antibody, EIF3-theta Antibody
Host	EIF3A antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Immunogen Species:	EIF3A / EIF3 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	EIF3A / EIF3 antibody was raised against a portion of amino acids 550-600 of human eIF3A was used as the immunogen.
Specificity:	Human EIF3A
Epitope:	aa550-600
Reactivity:	Human, Chimpanzee, Mouse, Rat, Bovine, Dog, Horse, Chicken
Purification:	Immunoaffinity purified
Presentation:	PBS, 0.05% BSA, 0.05% sodium azide
Recommended Storage:	Store at 4°C for short term applications. For long term storage, aliquot and store at -20°C.
Uses:	IHC - Paraffin (10 µg/ml), Western blot (0.25 - 1 µg/ml) (Optimal dilution to be determined by the researcher)
Size:	50 µg
Concentration:	0.5 mg/ml

Immunohistochemistry Image:



Anti-eIF3 antibody IHC of human kidney. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval. Antibody LS-B7957 dilution 10 ug/ml.

Western Blot Image:



Western blot of eIF3A in Daudi cell lysate in the 1) absence and 2) presence of immunizing peptide and 3) NIH 3T3 cell lysate using LS-B7957 at 0.5 ug/ml, 0.5 ug/ml and 0.25 ug/ml respectively.

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences