


Cytokeratin 5+6 Mouse anti-Human Monoclonal (D5/16 B4) Antibody - LS-B8554 - LSBio

CatalogID:	LS-B8554
Validation:	This antibody replaces catalog number LS-C87476. It has been validated for use in the following assays: IHC-P.
Target:	Keratin 5/Keratin 6
Host	Keratin 5/Keratin 6 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1,k
Clone Name:	D5/16 B4
Immunogen Species:	Cytokeratin 5+6 antibody was raised against Human
Antigen Type:	Purified protein
Immunogen:	Cytokeratin 5+6 antibody was raised against isolated cytokeatin 5.
Reactivity:	Human
Presentation:	0.09% sodium azide. Tissue Culture Supernatant
Recommended Storage:	+4°C, avoid freezing
Uses:	IHC - Paraffin (1:50) (Optimal dilution to be determined by the researcher)
Size:	250 µl
Concentration:	200 µg/ml

Immunohistochemistry Image:



Anti-Keratin 5/6 antibody IHC of human skin, epidermis. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval. Antibody LS-B8554 dilution 1:50.

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences