


HMGB1 / HMGB1 Mouse anti-Human Monoclonal (IMG19N15F4) Antibody - LS-B8674 - LSBio	
CatalogID:	LS-B8674
Validation:	This antibody replaces catalog number LS-C148831. It has been validated for use in the following assays: IHC-P.
Target:	high mobility group box 1 (HMGB1)
Synonyms:	HMGB1 Antibody, Amphoterin Antibody, High mobility group protein B1 Antibody, High mobility group box 1 Antibody, High mobility group protein 1 Antibody, HMG1 Antibody, HMG3 Antibody, SBP-1 Antibody, High-mobility group box 1 Antibody, HMG-1 Antibody
Host	HMGB1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1,k
Clone Name:	IMG19N15F4
Immunogen Species:	HMGB1 / HMGB1 antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	HMGB1 / HMGB1 antibody was raised against full-length recombinant human HMGB1 was used as the immunogen.
Specificity:	Human HMGB1
Reactivity:	Human, Mouse
Purification:	Protein G purified
Presentation:	PBS, 0.05% BSA, 0.05% sodium azide
Recommended Storage:	Store at 4°C for short term applications. For long term storage, aliquot and store at -20°C.
Usage Summary:	Flow: Intracellular.
Uses:	IHC - Paraffin (5 µg/ml), Western blot (1 - 3 µg/ml), Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	50 µg
Concentration:	0.5 mg/ml

Immunohistochemistry Image:



Anti-HMG1 / HMGB1 antibody IHC of human skin, epidermis. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval. Antibody LS-B8674 dilution 5 ug/ml.

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences