


Basigin / Emmprin / CD147 Mouse anti-Human Monoclonal (56.3_1D5) Antibody - LS-B8722 - LSBio	
CatalogID:	LS-B8722
Validation:	This antibody replaces catalog number LS-C155776. It has been validated for use in the following assays: IHC-P.
Target:	basigin (Ok blood group) (BSG)
Synonyms:	BSG Antibody, 5F7 Antibody, Basigin Antibody, Basigin (Ok blood group) Antibody, CD147 antigen Antibody, Collagenase stimulatory factor Antibody, EMMPRIN Antibody, OK Antibody, OK blood group antigen Antibody, TCSF Antibody, CD147 Antibody, M6 Antibody
Host	BSG antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2b
Clone Name:	56.3_1D5
Immunogen Species:	Basigin / Emmprin / CD147 antibody was raised against Human
Antigen Type:	Cells
Immunogen:	Basigin / Emmprin / CD147 antibody was raised against myeloid Leukemia cells (HL60)
Specificity:	Human Emmprin / CD147
Reactivity:	Human
Purification:	Protein A/G purified
Presentation:	PBS, 0.02% sodium azide
Recommended Storage:	Long term: -20°C; Short term: +4°C. Avoid repeat freeze-thaw cycles.
Uses:	IHC - Paraffin (10 µg/ml), ICC (10 µg/ml), Western blot (1 µg/ml) (Optimal dilution to be determined by the researcher)
Size:	50 µg

Immunohistochemistry Image:



Human Testis: Formalin-Fixed, Paraffin-Embedded (FFPE)

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences