



DCN / Decorin Rabbit anti-Human Polyclonal Antibody - LS-B8745 - LSBio	
<b>CatalogID:</b>	LS-B8745
<b>Validation:</b>	This antibody replaces catalog number LS-C146938. It has been validated for use in the following assays: IHC-P.
<b>Target:</b>	decorin (DCN)
<b>Synonyms:</b>	DCN Antibody, CSCD Antibody, Decorin proteoglycan Antibody, PG-S2 Antibody, PG40 Antibody, Proteoglycan core protein Antibody, SLRR1B Antibody, Small leucine-rich protein 1B Antibody, PGI Antibody, PGS2 Antibody, Bone proteoglycan II Antibody, Decorin Antibody, DSPG2 Antibody
<b>Host</b>	DCN antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Immunogen Species:</b>	DCN / Decorin antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	DCN / Decorin antibody was raised against synthetic peptide surrounding amino acid 130 of human Decorin
<b>Specificity:</b>	Human DCN / Decorin
<b>Reactivity:</b>	Human, Mouse, Rat, Pig, Rabbit
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 25% glycerol, 0.5% BSA, 0.01% Thimerosal
<b>Recommended Storage:</b>	Long term: -20°C; Short term: +4°C. Avoid repeat freeze-thaw cycles.
<b>Uses:</b>	IHC - Paraffin (5 µg/ml), Western blot (0.5 - 4 µg/ml) (Optimal dilution to be determined by the researcher)
<b>Size:</b>	50 µg
<b>Concentration:</b>	0.5 mg/ml

**Immunohistochemistry Image:**



Human Skin: Formalin-Fixed, Paraffin-Embedded (FFPE)

**Western Blot Image:**



Western blot of Decorin antibody in Jurkat cell lysate.

**Requested From:**

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences