


TPT1 / TCTP Mouse anti-Human Monoclonal Antibody - LS-B8953 - LSBio


CatalogID:	LS-B8953
Validation:	This antibody replaces catalog number LS-C137540. It has been validated for use in the following assays: IHC-P.
Target:	tumor protein, translationally-controlled 1 (TPT1)
Synonyms:	TPT1 Antibody, Histamine-releasing factor Antibody, Fortilin Antibody, p02 Antibody, p23 Antibody, HRF Antibody, TCTP Antibody
Host	TPT1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Immunogen Species:	TPT1 / TCTP antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	TPT1 / TCTP antibody was raised against recombinant protein corresponding to full length TPT1
Specificity:	Anti-TPT1 recognizes TPT1 expression in SW480 whole cell lysate.
Reactivity:	Human
Purification:	Affinity purified
Presentation:	Citrate-Tris-HCl buffer, pH 7.0, 0.02% Proclin 300
Recommended Storage:	Short term 4°C, long term aliquot and store at -20°C, avoid freeze thaw cycles.
Uses:	IHC - Paraffin (10 µg/ml), Western blot (1:1000), ELISA (Optimal dilution to be determined by the researcher)
Size:	50 µl
Concentration:	1 mg/ml

Immunohistochemistry Image:



Human Spleen: Formalin-Fixed, Paraffin-Embedded (FFPE)

Western Blot Image:



Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences