



| UMOD / Uromodulin Rabbit anti-Human Polyclonal (aa351-380) Antibody - LS-B9021 - LSBio |  |
|--|--|
| <b>CatalogID:</b>  | LS-B9021   |
| <b>Validation:</b>   | This antibody replaces catalog number LS-C160430. It has been validated for use in the following assays: IHC-P.  |
| <b>Target:</b>   | uromodulin (UMOD)  |
| <b>Synonyms:</b>   | UMOD Antibody, ADMCKD2 Antibody, FJHN Antibody, HNFJ Antibody, MCKD2 Antibody, THGP Antibody, THP Antibody, Uromucoid Antibody, Tamm-Horsfall glycoprotein Antibody, HNFJ1 Antibody, Uromodulin Antibody |
| <b>Host</b>  | UMOD antibody was produced in Rabbit   |
| <b>Clonality:</b>  | Polyclonal   |
| <b>Immunogen Species:</b>  | UMOD / Uromodulin antibody was raised against Human  |
| <b>Antigen Type:</b>   | Synthetic peptide  |
| <b>Immunogen:</b>  | UMOD / Uromodulin antibody was raised against kLH-conjugated synthetic peptide from internal region of human UMOD.   |
| <b>Specificity:</b>  | Human UMOD / Uromodulin  |
| <b>Epitope:</b>  | aa351-380  |
| <b>Reactivity:</b>   | Human  |
| <b>Purification:</b>   | Immunoaffinity purified  |
| <b>Presentation:</b>   | PBS, 0.09% sodium azide  |
| <b>Recommended Storage:</b>  | Maintain refrigerated at 2-8°C for up to 6 months. For long term storage store at -20°C.   |
| <b>Uses:</b>   | IHC - Paraffin (5 µg/ml), Western blot (1:1000) (Optimal dilution to be determined by the researcher)  |
| <b>Size:</b>   | 200 µl   |

**Immunohistochemistry Image:**



Human Kidney: Formalin-Fixed, Paraffin-Embedded (FFPE)

**Western Blot Image:**



UMOD Antibody (Center) Western blot of MDA-MB453 cell line lysates (35 ug/lane). This demonstrates the UMOD antibody detected the UMOD protein (arrow).

**Requested From:**

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences