



HOXA7 Rabbit anti-Rat Polyclonal (aa1-30) Antibody - LS-B9160 - LSBio	
CatalogID:	LS-B9160
Validation:	This antibody replaces catalog number LS-C162468. It has been validated for use in the following assays: IHC-P.
Target:	homeobox A7 (HOXA7)
Synonyms:	HOXA7 Antibody, Homeotic protein hox-a7 Antibody, Homeo box A7 Antibody, Homeobox A7 Antibody, Homeobox protein Hox 1.1 Antibody, Homeobox protein Hox-1A Antibody, HOX1 Antibody, ANTP Antibody, Homeobox protein Hox-A7 Antibody, HOX1.1 Antibody, HOX1A Antibody
Host	HOXA7 antibody was produced in Rabbit
Clonality:	Polyclonal
Immunogen Species:	HOXA7 antibody was raised against Rat
Antigen Type:	Synthetic peptide - KLH conjugated
Immunogen:	HOXA7 antibody was raised against KLH-conjugated synthetic peptide from N-terminal region of rat Hoxa7
Epitope:	aa1-30
Reactivity:	Rat, Human
Purification:	Immunoaffinity purified
Presentation:	PBS with 0.09% (w/v) sodium azide
Recommended Storage:	Maintain refrigerated at 2-8°C for up to 6 months. For long term storage store at -20°C.
Uses:	IHC - Paraffin (5 µg/ml), Western blot (1:1000) (Optimal dilution to be determined by the researcher)
Size:	200 µl

Immunohistochemistry Image:



Human, Skin: Formalin-Fixed Paraffin-Embedded (FFPE)

Western Blot Image:



Rat Hoxa7 Antibody (N-term) Western blot analysis in WiDr cell line lysates (35 ug/lane). This demonstrates the Hoxa7 antibody detected the Hoxa7 protein (arrow).

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences