


STAT3 Rabbit anti-Mouse Polyclonal (C-Terminus) Antibody - LS-B9434 - LSBio

CatalogID:	LS-B9434
Validation:	This antibody replaces catalog number LS-C51705. It has been validated for use in the following assays: IHC-P.
Target:	signal transducer and activator of transcription 3 (acute-phase response factor) (STAT3)
Synonyms:	STAT3 Antibody, Acute-phase response factor Antibody, APRF Antibody, DNA-binding protein APRF Antibody, HIES Antibody
Host	STAT3 antibody was produced in Rabbit
Clonality:	Polyclonal
Immunogen Species:	STAT3 antibody was raised against Mouse
Immunogen:	STAT3 antibody was raised against peptide mapping to the carboxy terminus of mouse Stat3 p92.
Specificity:	STAT3 p92 of mouse, rat, and human origin
Epitope:	C-Terminus
Reactivity:	Mouse, Human, Rat
Purification:	Immunoaffinity purified
Presentation:	PBS, 0.1% sodium azide, 0.2% gelatin.
Recommended Storage:	Store at +4°C. DO NOT FREEZE.
Uses:	IHC - Paraffin (5 µg/ml), Western blot (1:200), Immunoprecipitation (Optimal dilution to be determined by the researcher)
Size:	50 µg
Concentration:	0.2 mg/ml

Immunohistochemistry Image:



Human Skin: Formalin-Fixed, Paraffin-Embedded (FFPE)

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences