



DSC3 / Desmocollin 3 Rabbit anti-Human Polyclonal (aa821-850) Antibody - LS-B9474 - LSBio	
<b>CatalogID:</b>	LS-B9474
<b>Validation:</b>	This antibody replaces catalog number LS-C167540. It has been validated for use in the following assays: IHC-P.
<b>Target:</b>	desmocollin 3 (DSC3)
<b>Synonyms:</b>	DSC3 Antibody, CDHF3 Antibody, Cadherin family member 3 Antibody, Desmocollin-3 Antibody, Desmocollin 3 Antibody, DSC4 Antibody, HT-CP Antibody, Desmocollin-4 Antibody, DSC Antibody
<b>Family / Subfamily:</b>	Cadherin / not assigned-Cadherin
<b>Host</b>	DSC3 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Immunogen Species:</b>	DSC3 / Desmocollin 3 antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	DSC3 / Desmocollin 3 antibody was raised against kLH-conjugated synthetic peptide from C-terminal region of human DSC3.
<b>Specificity:</b>	Human DSC3
<b>Epitope:</b>	aa821-850
<b>Reactivity:</b>	Human, Mouse
<b>Purification:</b>	Immunoaffinity purified
<b>Presentation:</b>	PBS, 0.09% sodium azide
<b>Recommended Storage:</b>	Maintain refrigerated at 2-8°C for up to 6 months. For long term storage store at -20°C.
<b>Uses:</b>	IHC - Paraffin (1:100), Western blot (1:1000) (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl

**Immunohistochemistry Image:**



Human Placenta: Formalin-Fixed, Paraffin-Embedded (FFPE)

**Western Blot Image:**



DSC3 Antibody (C-term) Western blot of mouse spleen tissue lysates (35 ug/lane). This demonstrates the DSC3 antibody detected the DSC3 protein (arrow).

<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/24/2014	
© 2014 LifeSpan BioSciences	