


ACTN3 Rabbit anti-Human Polyclonal (aa409-437) Antibody - LS-B9666 - LSBio


CatalogID:	LS-B9666
Validation:	This antibody replaces catalog number LS-C156007. It has been validated for use in the following assays: IHC-P.
Target:	actinin, alpha 3 (ACTN3)
Synonyms:	ACTN3 Antibody, Actinin, alpha 3 Antibody, Alpha-actinin skeletal muscle Antibody, Alpha-actinin-3 Antibody, F-actin cross-linking protein Antibody
Host	ACTN3 antibody was produced in Rabbit
Clonality:	Polyclonal
Immunogen Species:	ACTN3 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	ACTN3 antibody was raised against kLH-conjugated synthetic peptide from internal region of human ACTN3.
Specificity:	Human ACTN3
Epitope:	aa409-437
Reactivity:	Human
Purification:	Immunoaffinity purified
Presentation:	PBS, 0.09% sodium azide
Recommended Storage:	Maintain refrigerated at 2-8°C for up to 6 months. For long term storage store at -20°C.
Uses:	IHC - Paraffin (1:100), Western blot (1:1000) (Optimal dilution to be determined by the researcher)
Size:	200 µl

Immunohistochemistry Image:



Human Skeletal Muscle: Formalin-Fixed, Paraffin-Embedded (FFPE)

Western Blot Image:



ACTN3 Antibody (Center) Western blot of MDA-MB231 cell line lysates (35 ug/lane). This demonstrates the ACTN3 antibody detected the ACTN3 protein (arrow).

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences