



CA3 / Carbonic Anhydrase III Rabbit anti-Human Polyclonal (N-Terminus) Antibody - LS-B9781 - LSBio	
CatalogID:	LS-B9781
Validation:	This antibody replaces catalog number LS-C97508. It has been validated for use in the following assays: IHC-P.
Target:	carbonic anhydrase III, muscle specific (CA3)
Synonyms:	CA3 Antibody, Carbonate dehydratase III Antibody, Carbonic anhydrase 3 Antibody, Carbonic anhydrase III Antibody, CA-III Antibody, CAIII Antibody, Car3 Antibody
Host	CA3 antibody was produced in Rabbit
Clonality:	Polyclonal
Immunogen Species:	CA3 / Carbonic Anhydrase III antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	CA3 / Carbonic Anhydrase III antibody was raised against kLH conjugated synthetic peptide selected from the N-terminal region of human CA3.
Epitope:	N-Terminus
Reactivity:	Human, Mouse
Purification:	Ammonium sulfate precipitation
Presentation:	PBS, 0.09% sodium azide
Recommended Storage:	Maintain refrigerated at 2-8°C for up to 6 months. For long term storage store at -20°C.
Uses:	IHC - Paraffin (1:50), Western blot (1:1000) (Optimal dilution to be determined by the researcher)
Size:	200 µl

Immunohistochemistry Image:




Human Skeletal Muscle: Formalin-Fixed, Paraffin-Embedded (FFPE)

Western Blot Image:



Western blot of CA3 antibody (N-term) in mouse kidney tissue lysates (35 ug/lane). CA3 (arrow) was detected using the purified Pab.

Western Blot Image:



Western blot of CA3 antibody (N-term) in HepG2 cell line lysates (35 ug/lane). CA3 (arrow) was detected using the purified Pab.

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences