



ALK3 / BMPR1A Rabbit anti-Human Polyclonal (aa169-197) Antibody - LS-B9804 - LSBio	
CatalogID:	LS-B9804
Validation:	This antibody replaces catalog number LS-C164304. It has been validated for use in the following assays: IHC-P.
Target:	bone morphogenetic protein receptor, type IA (BMPR1A)
Synonyms:	BMPR1A Antibody, 10q23del Antibody, ALK3 Antibody, ALK-3 Antibody, BMP type -1A receptor Antibody, BMPR-1A Antibody, CD292 Antibody, CD292 antigen Antibody, Activin receptor-like kinase 3 Antibody, ACVRLK3 Antibody, SKR5 Antibody
Family / Subfamily:	Protein Kinase / Activin/TGFb
Host	BMPR1A antibody was produced in Rabbit
Clonality:	Polyclonal
Immunogen Species:	ALK3 / BMPR1A antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	ALK3 / BMPR1A antibody was raised against kLH-conjugated synthetic peptide from internal region of human Mouse Bmpr1a.
Specificity:	Human ALK3 / BMPR1A
Epitope:	aa169-197
Reactivity:	Human, Mouse
Purification:	Immunoaffinity purified
Presentation:	PBS, 0.09% sodium azide
Recommended Storage:	Maintain refrigerated at 2-8°C for up to 6 months. For long term storage store at -20°C.
Uses:	IHC - Paraffin (1:100), Western blot (1:1000) (Optimal dilution to be determined by the researcher)
Size:	200 µl

Immunohistochemistry Image:



Human Small Intestine: Formalin-Fixed, Paraffin-Embedded (FFPE)

Western Blot Image:



Mouse Bmpr1a Antibody (Center) Western blot of SK-BR-3 cell line lysates (35 ug/lane). This demonstrates the Bmpr1a antibody detected the Bmpr1a protein (arrow).

Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/24/2014	
© 2014 LifeSpan BioSciences	