


TMX3 Rabbit anti-Human Polyclonal (aa402-431) Antibody - LS-B9820 - LSBio


CatalogID:	LS-B9820
Validation:	This antibody replaces catalog number LS-C156248. It has been validated for use in the following assays: IHC-P.
Target:	thioredoxin-related transmembrane protein 3 (TMX3)
Synonyms:	TMX3 Antibody, KIAA1830 Antibody, TXNDC10 Antibody, PDIA13 Antibody
Host	TMX3 antibody was produced in Rabbit
Clonality:	Polyclonal
Immunogen Species:	TMX3 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	TMX3 antibody was raised against kLH-conjugated synthetic peptide from C-terminal region of human TMX3.
Specificity:	Human TMX3
Epitope:	aa402-431
Reactivity:	Human
Purification:	Immunoaffinity purified
Presentation:	PBS, 0.09% sodium azide
Recommended Storage:	Maintain refrigerated at 2-8°C for up to 6 months. For long term storage store at -20°C.
Uses:	IHC - Paraffin (1:100), Western blot (1:1000) (Optimal dilution to be determined by the researcher)
Size:	200 µl

Immunohistochemistry Image:



Human Liver: Formalin-Fixed, Paraffin-Embedded (FFPE)

Western Blot Image:



TMX3 Antibody (C-term) Western blot of A2058 cell line lysates (35 ug/lane). This demonstrates the TMX3 antibody detected the TMX3 protein (arrow).

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences