


MFAP4 Rabbit anti-Human Polyclonal (aa225-255) Antibody - LS-B9892 - LSBio

CatalogID:	LS-B9892
Validation:	This antibody replaces catalog number LS-C164160. It has been validated for use in the following assays: IHC-P.
Target:	microfibrillar-associated protein 4 (MFAP4)
Synonyms:	MFAP4 Antibody
Host	MFAP4 antibody was produced in Rabbit
Clonality:	Polyclonal
Immunogen Species:	MFAP4 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	MFAP4 antibody was raised against kLH-conjugated synthetic peptide from C-terminal region of human MFAP4.
Specificity:	Human MFAP4
Epitope:	aa225-255
Reactivity:	Human, Mouse
Purification:	Immunoaffinity purified
Presentation:	PBS, 0.09% sodium azide
Recommended Storage:	Maintain refrigerated at 2-8°C for up to 6 months. For long term storage store at -20°C.
Uses:	IHC - Paraffin (1:50 - 1:100), Western blot (1:1000), Flow Cytometry (1:10 - 1:50) (Optimal dilution to be determined by the researcher)
Size:	200 µl


Immunohistochemistry Image:



Human Tonsil, Collagenous Stroma: Formalin-Fixed, Paraffin-Embedded (FFPE)


Flow Cytometry Image:

MCF-7



MFAP4 Antibody (C-term) flow cytometry of MCF-7 cells (right histogram) compared to a negative control cell (left histogram). FITC-conjugated goat-anti-rabbit secondary antibodies were used for the analysis.

Western Blot Image:



Western blot of MFAP4 Antibody (C-term) in mouse lung tissue lysates (35 ug/lane). MFAP4 (arrow) was detected using the purified Pab.

Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/24/2014	
© 2014 LifeSpan BioSciences	