

E. coli O157:H7 Mouse anti-E. coli Monoclonal (FITC) Antibody - LS-C103125 - LSBio

CatalogID:	LS-C103125
Target:	E. coli O157:H7
Host	E. coli O157:H7 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2b
Conjugations:	Fluorescein (FITC)
Immunogen Species:	E. coli O157:H7 antibody was raised against E. coli
Specificity:	O157
Reactivity:	E. coli
Purification:	Protein A purified
Presentation:	0.01 M PBS, pH 7.2, 0.1% sodium azide, bovine serum albumin at 10 mg/ml
Recommended Storage:	Store at -20°C. Aliquot to avoid freeze/thaw cycles.
Uses:	Immunofluorescence (1:10 - 1:50), ELISA (1:20 - 1:200) (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences