

**EBV / Epstein-Barr Virus Mouse anti-Epstein-Barr Virus Monoclonal (FITC) Antibody - LS-C103127
 - LSBio**

CatalogID:	LS-C103127
Target:	EBV / Epstein-Barr Virus
Host	Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Epstein-Barr Virus
Specificity:	gp 250/350
Reactivity:	Epstein-Barr Virus
Purification:	Protein A purified
Presentation:	0.01 M PBS, pH 7.2, 0.1% sodium azide, bovine serum albumin at 10 mg/ml
Recommended Storage:	Store at -20°C. Aliquot to avoid freeze/thaw cycles.
Uses:	Immunofluorescence (1:10 - 1:50), Western blot (1:10 - 1:50) (Optimal dilution to be determined by the researcher)
Size:	100 µg
Concentration:	0.1 mg/ml
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences