

Giardia lamblia Mouse anti-Giardia lamblia Monoclonal (FITC) Antibody - LS-C103128 - LSBio

CatalogID:	LS-C103128
Target:	Giardia lamblia
Host	Giardia lamblia antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG3
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Giardia lamblia antibody was raised against Giardia lamblia
Specificity:	cysts
Reactivity:	Giardia lamblia
Purification:	Protein A purified
Presentation:	0.01 M PBS, pH 7.2, 0.1% sodium azide, bovine serum albumin at 10 mg/ml
Recommended Storage:	Store at -20°C. Aliquot to avoid freeze/thaw cycles.
Uses:	Immunofluorescence (1:10 - 1:50), ELISA (1:20 - 1:200) (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences