

Influenza A Virus Nucleoprotein Mouse anti-Influenza Virus Monoclonal (FITC) Antibody - LS-C103137 - LSBio	
<b>CatalogID:</b>	LS-C103137
<b>Target:</b>	Influenza A Virus Nucleoprotein
<b>Host</b>	Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	Influenza Virus
<b>Specificity:</b>	NP
<b>Reactivity:</b>	Influenza Virus
<b>Purification:</b>	Protein A purified
<b>Presentation:</b>	0.01 M PBS, pH 7.2, 0.1% sodium azide, bovine serum albumin at 10 mg/ml
<b>Recommended Storage:</b>	Store at -20°C. Aliquot to avoid freeze/thaw cycles.
<b>Uses:</b>	Immunofluorescence (1:10 - 1:50), ELISA (1:20 - 1:200) (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan
<b>Laboratory Reagent For In Vitro Research Use Only</b> <b>Not for resale without prior written consent from LifeSpan BioSciences, Inc.</b> Created on 9/24/2014 © 2014 LifeSpan BioSciences	