

Candida albicans Rabbit anti-Candida albicans Polyclonal (FITC) Antibody - LS-C103355 - LSBio	
CatalogID:	LS-C103355
Target:	Candida albicans
Host	Rabbit
Clonality:	Polyclonal
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Candida albicans
Antigen Type:	Cells
Immunogen:	Candida albicans, type A.
Specificity:	Recognizes numerous proteins in a soluble C. albicans extract (IEP). This antibody does cross react with other yeasts.
Reactivity:	Candida albicans
Purification:	Purified IgG
Presentation:	0.01 M PBS, pH 7.2, 0.1% sodium azide
Recommended Storage:	Long term: -20°C; Short term: +4°C; Avoid freeze-thaw cycles.
Uses:	Immunofluorescence (1:1200) (Optimal dilution to be determined by the researcher)
Size:	1 ml
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 9/24/2014 © 2014 LifeSpan BioSciences	