

**HSV / Herpes Simplex Virus Rabbit anti-Herpes Simplex Virus Polyclonal (FITC) Antibody - LS-C103373 - LSBio**

<b>CatalogID:</b>	LS-C103373
<b>Target:</b>	Herpes Simplex Virus (HSV-1/HSV-2)
<b>Host</b>	Rabbit
<b>Clonality:</b>	Polyclonal
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	Herpes Simplex Virus
<b>Specificity:</b>	types 1 &2
<b>Reactivity:</b>	Herpes Simplex Virus
<b>Purification:</b>	Purified IgG
<b>Presentation:</b>	0.01 M PBS, pH 7.2, 0.1% sodium azide
<b>Recommended Storage:</b>	Long term: -20°C; Short term: +4°C; Avoid freeze-thaw cycles.
<b>Uses:</b>	Immunofluorescence (1:2000) (Optimal dilution to be determined by the researcher)
<b>Size:</b>	1 ml
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences