

HSV-2 / Herpes Simplex Virus 2 Sheep anti-Herpes Simplex Virus Polyclonal (FITC) Antibody - LS-C103374 - LSBio	
CatalogID:	LS-C103374
Target:	Herpes Simplex Virus 2 (HSV-2)
Host	Sheep
Clonality:	Polyclonal
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Herpes Simplex Virus
Specificity:	type 2
Reactivity:	Herpes Simplex Virus
Purification:	Purified IgG
Presentation:	0.01 M PBS, pH 7.2, 0.1% sodium azide
Recommended Storage:	Long term: -20°C; Short term: +4°C; Avoid freeze-thaw cycles.
Uses:	Immunofluorescence (1:4000) (Optimal dilution to be determined by the researcher)
Size:	1 ml
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/24/2014	
© 2014 LifeSpan BioSciences	