

Mycoplasma pneumoniae Rabbit anti-Mycoplasma pneumoniae Polyclonal (FITC) Antibody - LS-C103389 - LSBio

CatalogID:	LS-C103389
Target:	Mycoplasma pneumoniae
Host	Rabbit
Clonality:	Polyclonal
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Mycoplasma pneumoniae
Antigen Type:	Cells
Immunogen:	Purified M. pneumoniae.
Specificity:	all antigens
Reactivity:	Mycoplasma pneumoniae
Purification:	Purified IgG
Presentation:	0.01 M PBS, pH 7.2, 0.1% sodium azide
Recommended Storage:	Long term: -20°C; Short term: +4°C; Avoid freeze-thaw cycles.
Uses:	Immunofluorescence (1:1000) (Optimal dilution to be determined by the researcher)
Size:	1 ml
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences