


CD45.1 Mouse anti-Mouse Monoclonal (FITC) (A20) Antibody - LS-C105864 - LSBio


CatalogID:	LS-C105864
Target:	CD45.1
Host	CD45.1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a,k
Clone Name:	A20
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD45.1 antibody was raised against Mouse
Immunogen:	CD45.1 antibody was raised against mouse CD45.1
Reactivity:	Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 150 mM sodium chloride, 0.09% sodium azide
Recommended Storage:	Store at +4°C. Do not freeze. Product is photosensitive and should be protected from light.
Usage Summary:	The A20 antibody has been tested by flow cytometric analysis of mouse splenocyte suspensions. This can be used at less than or equal to 0.5 ug per test. A test is defined as the amount (ug) of antibody that will stain a cell sample in a final volume of 100 ul. Cell number should be determined empirically but can range from 10 ⁵ to 10 ⁸ cells/test. It is recommended that the antibody be carefully titrated for optimal performance in the assay of interest.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	50 µg or 100 µg or 500 µg

Flow Cytometry Image:



Staining of BALB/c splenocytes with 0.25 ug FITC Mouse IgG isotype control or 0.25 ug FITC anti-mouse CD45.1 (A20) (colored histogram). Total viable cells were used for analysis.

Flow Cytometry Image:



Staining of SJL splenocytes with 0.25 ug FITC Mouse IgG isotype control or 0.25 ug FITC anti-mouse CD45.1 (A20) (colored histogram). Total viable cells were used for analysis.

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences