


CD4 Rat anti-Mouse Monoclonal (FITC) (RM4-4) Antibody - LS-C106097 - LSBio

CatalogID:	LS-C106097
Target:	CD4 molecule
Synonyms:	CD4 Antibody, CD4 antigen Antibody, CD4 receptor Antibody, CD4mut Antibody, CD4 antigen (p55) Antibody, CD4 molecule Antibody
Host	CD4 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2b,k
Clone Name:	RM4-4
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD4 antibody was raised against Mouse
Immunogen:	CD4 antibody was raised against mouse CD4
Reactivity:	Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 150 mM sodium chloride, 0.09% sodium azide
Recommended Storage:	Store at +4°C. Do not freeze. Product is photosensitive and should be protected from light.
Usage Summary:	This RM4-4 antibody has been tested by flow cytometric analysis of mouse splenocytes. This can be used at less than or equal to 0.125 ug per test. A test is defined as the amount (ug) of antibody that will stain a cell sample in a final volume of 100 ul. Cell number should be determined empirically but can range from 10 ⁵ to 10 ⁸ cells/test. It is recommended that the antibody be carefully titrated for optimal performance in the assay of interest.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	50 µg or 100 µg or 500 µg

Flow Cytometry Image:



Staining of C57Bl/6 splenocytes with .125 ug of FITC Rat IgG2b Isotype Control (open histogram) or .06 ug of FITC anti-mouse CD4 (RM4-4) (colored histogram). Cells in the lymphocyte gate were used for analysis.

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences