

TLR3 Mouse anti-Human Monoclonal (PE) Antibody - LS-C10692 - LSBio

CatalogID:	LS-C10692
Target:	toll-like receptor 3 (TLR3)
Synonyms:	TLR3 Antibody, CD283 Antibody, Toll-like receptor 3 Antibody, CD283 antigen Antibody, IIAE2 Antibody
Family / Subfamily:	Toll-like Receptor / not assigned-Toll-like Receptor
Host	TLR3 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	TLR3 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	TLR3 antibody was raised against kLH-conjugated synthetic peptide of human TLR3.
Specificity:	Recognizes human TLR3 (Toll Like Receptor 3).
Reactivity:	Human
Purification:	Protein G purified
Presentation:	PBS, 0.02% sodium azide.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Suitable for use in Western Blot and Flow Cytometry. Western Blot: 1-3 ug/ml. Flow Cytometry (Intercellular): 2-5 ug/1x1E6 cells. Positive control: Human Intestine or Ramos.
Uses:	Western blot (1 - 3 µg/ml), Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences