

TLR4 Mouse anti-Human Monoclonal (PE) Antibody - LS-C10729 - LSBio	
<b>CatalogID:</b>	LS-C10729
<b>Target:</b>	toll-like receptor 4 (TLR4)
<b>Synonyms:</b>	TLR4 Antibody, CD284 Antibody, CD284 antigen Antibody, Homolog of Drosophila toll Antibody, TLR-4 Antibody, TOLL Antibody, ARMD10 Antibody, Toll-like receptor 4 Antibody, HToll Antibody
<b>Family / Subfamily:</b>	Toll-like Receptor / not assigned-Toll-like Receptor
<b>Host</b>	TLR4 antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Immunogen Species:</b>	TLR4 antibody was raised against Human
<b>Antigen Type:</b>	Cells
<b>Immunogen:</b>	TLR4 antibody was raised against ba/F3 cell line expressing TLR4.
<b>Specificity:</b>	Recognizes the human Toll like receptor 4 (TLR4) cell surface antigen. Species cross-reactivity: Monkey, Rhesus monkey.
<b>Reactivity:</b>	Human, Monkey
<b>Purification:</b>	Protein G purified
<b>Presentation:</b>	PBS, 0.09% sodium azide, 1% BSA, 40% glycerol.
<b>Recommended Storage:</b>	+4°C, avoid freezing
<b>Usage Summary:</b>	Suitable for use in Flow Cytometry. Flow Cytometry: Neat. Use 10 ul of the suggested working dilution to label 1E6 cells or 100 ul whole blood.
<b>Uses:</b>	Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	50 tst
<b>Requested From:</b>	Japan
<b>Laboratory Reagent For In Vitro Research Use Only</b>	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/24/2014	
© 2014 LifeSpan BioSciences	