


CD16 + CD32 Rat anti-Mouse Monoclonal (PE, Cy7) (93) Antibody - LS-C107347 - LSBio

CatalogID:	LS-C107347
Target:	CD16/CD32 (FcγR3/FcγR2)
Host	CD16/CD32 (FcγR3/FcγR2) antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2a,I
Clone Name:	93
Conjugations:	Phycoerythrin (PE), Cy7
Immunogen Species:	CD16 + CD32 antibody was raised against Mouse
Immunogen:	CD16 + CD32 antibody was raised against mouse CD16/CD32 (FcγR3/FcγR2)
Reactivity:	Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, ≤0.09% sodium azide, protein stabilizer
Recommended Storage:	Store at +4°C. Do not freeze. Product is photosensitive and should be protected from light.
Usage Summary:	The 93 antibody has been tested by flow cytometric analysis of mouse splenocyte suspensions. This can be used at less than or equal to 0.25 ug per test. A test is defined as the amount (ug) of antibody that will stain a cell sample in a final volume of 100 ul. Cell number should be determined empirically but can range from 10 ⁵ to 10 ⁸ cells/test. It is recommended that the antibody be carefully titrated for optimal performance in the assay of interest.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	50 µg or 100 µg

Flow Cytometry Image:



Staining of BALB/c splenocytes with staining buffer (autofluorescence) (open histogram) or 0.125 ug PE-Cy7 anti-mouse CD16/32 (93) (colored histogram). Total cells were used for analysis.

Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/24/2014	
© 2014 LifeSpan BioSciences	