


SOX2 Mouse anti-Human Monoclonal (10F10C9) Antibody - LS-C108882 - LSBio	
CatalogID:	LS-C108882
Target:	SRY (sex determining region Y)-box 2 (SOX2)
Synonyms:	SOX2 Antibody, MCOPS3 Antibody, Transcription factor SOX2 Antibody, ANOP3 Antibody, Transcription factor SOX-2 Antibody, SRY-related HMG-box gene 2 Antibody
Family / Subfamily:	Transcription factor
Host	SOX2 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	10F10C9
Immunogen Species:	SOX2 antibody was raised against Human
Immunogen:	SOX2 antibody was raised against purified recombinant fragment of SOX2(aa1-170) expressed in E. coli.
Specificity:	Human SOX2
Reactivity:	Human
Purification:	Ascites
Presentation:	Ascites fluid, 0.03% sodium azide
Recommended Storage:	Long term: -20°C; Short term: +4°C. Avoid repeat freeze-thaw cycles.
Uses:	Western blot (1:500 - 1:2000), ELISA (1:10000) (Optimal dilution to be determined by the researcher)
Uses Not Recommended:	IHC - Paraffin
Size:	100 µl

Western Blot Image:



Western blot of SOX2 monoclonal antibody against truncated Trx-SOX2 recombinant protein(1) and truncated MBP-SOX2(aa1-170) recombinant protein(2).

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences