


POU5F1 / OCT4 Mouse anti-Human Monoclonal (1D10H6) Antibody - LS-C109108 - LSBio	
CatalogID:	LS-C109108
Target:	POU class 5 homeobox 1 (POU5F1)
Synonyms:	POU5F1 Antibody, Oct-4 Antibody, OCT3 Antibody, OCT4 Antibody, OTF4 Antibody, POU class 5 homeobox 1 Antibody, Octamer-binding protein 3 Antibody, Octamer-binding protein 4 Antibody, OTF-3 Antibody, Oct-3 Antibody, OTF3 Antibody
Host	POU5F1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	1D10H6
Immunogen Species:	POU5F1 / OCT4 antibody was raised against Human
Immunogen:	POU5F1 / OCT4 antibody was raised against purified recombinant fragment of Oct4(aa193-360) expressed in E. coli.
Specificity:	Human POU5F1
Reactivity:	Human
Purification:	Ascites
Presentation:	Ascites fluid, 0.03% sodium azide
Recommended Storage:	Long term: -20°C; Short term: +4°C. Avoid repeat freeze-thaw cycles.
Uses:	Western blot (1:500 - 1:2000), ELISA (Optimal dilution to be determined by the researcher)
Uses Not Recommended:	IHC - Paraffin
Size:	100 µl

Western Blot Image:



Western blot of anti-Oct4 monoclonal antibody against recombinant Oct4 protein with Trx tag(1).

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences