

SELPLG / PSGL-1 / CD162 Mouse anti-Human Monoclonal (FITC) (TC2) Antibody - LS-C112524 - LSBio	
CatalogID:	LS-C112524
Target:	selectin P ligand (SELPLG)
Synonyms:	SELPLG Antibody, CD162 antigen Antibody, CLA Antibody, CD162 Antibody, PSGL1 Antibody, PSGL-1 Antibody, Selectin P ligand Antibody
Host	SELPLG antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	TC2
Conjugations:	Fluorescein (FITC)
Immunogen Species:	SELPLG / PSGL-1 / CD162 antibody was raised against Human
Specificity:	The antibody TC2 reacts with CD162, a 220 kD type I integral membrane protein expressed as disulfide-linked homodimer (sialomucin family). CD162 is present on the most peripheral blood T lymphocytes, monocytes, granulocytes; it is also expressed on a subpopulation of B lymphocytes and CD34+ bone marrow cells.
Reactivity:	Human
Purification:	Protein A purified
Presentation:	15 mM sodium azide, 0.2% BSA, PBS
Recommended Storage:	Store at +4°C. Do not freeze. Product is photosensitive and should be protected from light.
Usage Summary:	The reagent is designed for Flow Cytometry analysis of human blood cells using 20 ul reagent / 100 ul of whole blood or 10 ⁶ cells in a suspension. The content of a vial (2 ml) is sufficient for 100 tests.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/24/2014	
© 2014 LifeSpan BioSciences	