

BCL2L1 / BCL-XL Mouse anti-Human Monoclonal (PE) Antibody - LS-C12128 - LSBio

CatalogID:	LS-C12128
Target:	BCL2-like 1 (BCL2L1)
Synonyms:	BCL2L1 Antibody, Apoptosis regulator Bcl-X Antibody, Bcl-xL Antibody, BCL-XL/S Antibody, BCLXS Antibody, Bcl-2-like protein 1 Antibody, BCL2-like 1 Antibody, BCL2L Antibody, Bcl-X Antibody, Bcl-xS Antibody, BCLXL Antibody, PPP1R52 Antibody, Bcl2-L-1 Antibody, BCLX Antibody
Family / Subfamily:	Apoptosis
Host	BCL2L1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG3
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	BCL2L1 / BCL-XL antibody was raised against Human
Immunogen:	BCL2L1 / BCL-XL antibody was raised against recombinant Bcl-xS.
Specificity:	Recognizes human Bcl-xL (Mr 29kD).
Reactivity:	Human
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Immunoprecipitation. Immunohistochemistry. Flow cytometry. Western Blot. Flow Cytometry: 0.3 ug/1E6 cells.
Uses:	IHC, Western blot, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences