

**BrdU Mouse Monoclonal Antibody - LS-C121437 - LSBio**

<b>CatalogID:</b>	LS-C121437
<b>Target:</b>	Bromodeoxyuridine (BrdU)
<b>Host</b>	Bromodeoxyuridine (BrdU) antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG1
<b>Immunogen:</b>	BrdU antibody was raised against bromodeoxyuridine (BSA).
<b>Specificity:</b>	Recognizes the thymidine analog bromodeoxyuridine (BrdU).
<b>Purification:</b>	Protein G purified
<b>Presentation:</b>	PBS, pH 7.4, 0.09% sodium azide
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. For long-term storage and to avoid repeated freezing and thawing, aliquot and store at -20°C. Aliquots are stable for at least 12 months at -20°C
<b>Usage Summary:</b>	Suitable for detecting incorporated BrdU in a wide variety of cell types and is suitable for use on tissue sections in double-labeling techniques. Applications: Suitable for use in Flow Cytometry and Immunohistochemistry (frozen, paraffin). Flow Cytometry: 1:25-1:100; 10 ul labels 10 <sup>6</sup> cells in 100ul.
<b>Uses:</b>	IHC - Paraffin, IHC - Frozen, Flow Cytometry (1:25 - 1:100) (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/24/2014

© 2014 LifeSpan BioSciences